

# PRODUCT BROCHURE



# Mintti Breeze

---

## Ear Thermometer

Mintti Health Australia specializes in advanced health monitoring solutions, enhancing wellness through advanced medical technology in Australia and New Zealand.

☎ 1300 175 337

✉ [info@minttihealthaustralia.com.au](mailto:info@minttihealthaustralia.com.au)

🌐 [www.minttihealthaustralia.com.au](http://www.minttihealthaustralia.com.au)

# Ear Thermometer - Mintti Breeze

---

## Introduction:

Mintti Breeze infrared forehead and ear thermometer integrates with the latest built-in smart chip and AI algorithm. It can measure temperature with the highest degree of medical accuracy for all ages. Whether you are at home or work as a physician, Mintti Breeze allows you to monitor your temperature easily.

To meet the users' different needs and circumstances, Mintti Breeze can measure temperatures with contact or contact-free methods. The contact-free method with the silent mode is ideal for infants without waking them up while taking the temperature. You can also easily switch to ear reading to ensure more accurate reading on older children and adults.

With the world's first front-facing LED display, you don't have to turn Mintti Breeze over to check the reading as the display is right in front of you. It can make taking temperature so convenient. The need for self-measuring becomes more comfortable and much easier.

# Ear Thermometer - Mintti Breeze

---

## Specification:

<b>Product Name</b>	Mintti Breeze Infrared Thermometer
<b>Product Size</b>	161.5mm *35.8mm *55.9mm(length x width x height)
<b>Product Weight</b>	62 g (without battery)
<b>Measurement Temperature</b>	Ear Mode 34.0°C~42.4°C (93.2°F~108.3°F) Forehead Mode 34.0°C~42.4°C (93.2°F~108.3°F)
<b>Measurement Accuracy</b>	Ear Mode $\pm 0.2^{\circ}\text{C}/0.4^{\circ}\text{F}$ Forehead Mode $\pm 0.2^{\circ}\text{C}/0.4^{\circ}\text{F}$
<b>Bluetooth</b>	4.0 Enabled, Smart APP
<b>Memory Storage</b>	Up to 32 Readings
<b>Storage Tem. &amp; Hum</b>	-20° ~ 55°C(-4°F to 131°F) & $\leq 90\%$ moisture noncondensing
<b>Power Supply</b>	DC1.5V X 2AAA
<b>Power Consumption</b>	When off $\leq 10$ UW When measuring $\leq 200$ MW
<b>Auto Power-off</b>	30 seconds
<b>Product Accessories</b>	AAA batteries x2 and User Manual x1 Maintenance Service Card x 1 and QC Certificate x 1

# Ear Thermometer - Mintti Breeze

---

## Features:

### 01: **Fast & Accurate Measuring**

Combining cutting-edge sensors with a smart chip and optimized AI algorithms, Mintti Breeze can prevent false readings and offer the highest degree of accuracy within one second.

---

### 02: **Intelligent Performance & User Friendly**

It features a built-in fever alarm to prevent guesswork and an auto shutdown for power saving. With one-touch operation, it is simple to use even for children and elderly people.

---

### 03: **Easy-to-Read Front Facing LED Display**

With the world's first front-facing display, it makes self-measuring so easy. There is no need to turn it over to check the temperature as it is right in front of you.

---

### 04: **More Durable with Long-lasting Battery**

The long-lasting battery with approximately 3000 readings per year and 200 minutes of constant use makes it more reliable and durable to use.

---

### 05: **Preventing Infection with the Contact-free Forehead Measuring Mode**

The contact-free forehead measuring mode can prevent infection and cross-contamination, minimizing the spread of disease among use

---

### 06: **Ergonomic and Elegant Patented Design**

The ergonomic design of Mintti Breeze makes taking temperature more comfortable and much easier. Its unique upright patented design makes it look more like artwork rather than a thermometer.

---